



T:SLIM X2 INSULIN PUMP


# New Features for Software Version 7.7

## Pump Sounds

The pump will keep its existing tones for Pump Alarms, Malfunctions, and Alerts. The notification for Reminders now has a new tone.




Priority	Sound	Event
High	3 sequences of 3 beeps	Pump Alarms and Malfunctions
Medium	2 sequences of 3 beeps	Pump Alerts that impact insulin delivery
Low	1 sequence of 3 beeps	Pump Alerts that have no insulin delivery impact
<b>(NEW)</b> Informational	2 short beeps	Informal notifications (Reminders)



## Sensor Expiring Soon

Your pump tells you how much time you have left until your continuous glucose monitoring (CGM, sold separately) sensor session is complete. The Sensor Expiring Soon screen will now show at 24 hours remaining instead of at six hours remaining. It will also show at two hours remaining and 30 minutes remaining before your session ends.



## Auto-Off Alarm

Your pump can stop insulin delivery and alert you if there has been no interaction with the pump within a specified period of time. The default for this alarm is now preset to off.


Updating your current pump will not change your existing setting and you can still set it anywhere between 5-24 hours, or off.

If you are not using a CGM, we recommend keeping this Alert set to on.

### Responsible Use of Control-IQ Technology


Even with advanced systems such as the t:slim X2 insulin pump with Control-IQ technology, you are still responsible for actively managing your diabetes. Control-IQ technology does not prevent all high and low blood glucose events. The system is designed to help reduce glucose variability, but it requires your accurate input of information, such as meals and periods of sleep or exercise. Control-IQ technology will not function as intended unless you use all system components, including your CGM, infusion sets and pump cartridges, as instructed. Importantly, the system cannot adjust your insulin dosing if the pump is not receiving CGM readings. Because there are situations and emergencies that the system may not be capable of identifying or addressing, always pay attention to your symptoms and treat according to your healthcare provider's recommendations.

## Control-IQ technology updates



### Exercise Timer

You can now set a timer for a particular exercise duration so you do not have to remember when to stop the Activity when you are done exercising. You can set the timer between 30 minutes and eight hours.



### Activity Transition

You no longer need to stop an active Activity before switching to a different Activity. For example, if Sleep is active, you can start Exercise simply by tapping Exercise. You do not need to manually stop the Sleep Activity before tapping Exercise.

**Important Safety Information:** The t:slim X2 insulin pump with Control-IQ technology (the System) consists of the t:slim X2 insulin pump, which contains Control-IQ technology, and a compatible continuous glucose monitor (CGM, sold separately). The t:slim X2 insulin pump is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in people requiring insulin. The t:slim X2 insulin pump can be used solely for continuous insulin delivery and as part of the System. When used with a compatible CGM, the System can be used to automatically increase, decrease, and suspend delivery of basal insulin based on CGM sensor readings and predicted glucose values. The System can also deliver correction boluses when the glucose value is predicted to exceed a predefined threshold. The pump and the System are indicated for use in individuals six years of age and greater. The pump and the System are intended for single user use. The pump and the System are indicated for use with U-100 insulin only.

**WARNING:** Control-IQ technology should not be used by anyone under the age of six years old. It should also not be used in users who require less than 10 units of insulin per day or who weigh less than 25 kilograms.

The System is not indicated for use in pregnant women, people on dialysis, or critically ill users. Do not use the System if using hydroxyurea.

Users of the pump and the System must: be willing and able to use the insulin pump, CGM, and all other system components in accordance with their respective instructions for use; test blood glucose levels as recommended by their healthcare provider; demonstrate adequate carb-counting skills; maintain sufficient diabetes self-care skills; see healthcare provider(s) regularly; and have adequate vision and/or hearing to recognize all functions of the pump, including alerts, alarms, and reminders. The t:slim X2 pump and the CGM transmitter and sensor must be removed before MRI, CT, or diathermy treatment. Visit [tandemdiabetes.com/safetyinfo](https://tandemdiabetes.com/safetyinfo) for additional important safety information.

© 2023 Tandem Diabetes Care, Inc. All rights reserved. Tandem Diabetes Care, the Tandem logo, Control-IQ, and t:slim X2 are either registered trademarks or trademarks of Tandem Diabetes Care, Inc. in the United States and/or other countries. ML-1012180\_A






#### Customer Service

UK 0800 0121 560  
IRE 1800 12 4912

t:simulator App  
A free virtual pump demo



 [makingdiabeteseasier.com/uk](https://makingdiabeteseasier.com/uk)  
 @ALDiabetesUK  
 @makingdiabeteseasier\_uk

Air Liquide Healthcare is a world leader in medical gases, home healthcare, hygiene products and healthcare specialty ingredients. It aims to provide customers in the continuum of care from hospital to home with medical products, specialty ingredients and services that contribute to protecting vulnerable lives.